

Human Endothelin 1 (ET-1) ELISA Kit

Cat #: orb1807555 (manual)

For research use only. Not intended for diagnostic use.

Product Features

Intended Use

This ELISA kit applies to the in vitro quantitative determination of Human ET-1 concentrations in serum, plasma and other biological fluids.

Sensitivity	0.75 pg/mL
Detection Range	1.25-80pg/mL
Specificity	This kit recognizes Human ET-1 in samples. No significant cross-reactivity or interference between Human ET-1 and analogs was observed.
Repeatability	Coefficient of variation is $\leq 10\%$

Kit Components

3. Samples should be assayed within 7 days when stored at 2-8°C, otherwise samples must be divided up and stored at -20°C (<

Assay Procedure Summary

Calculation of Results

rate is divided into spiked recovery and

||
||
||

||
||
||

||
||
||